Industrial  
evaporation

## Betriebsanleitung Operating instructions Notice d'instructions

Hei-VAP Industrial

 **heidolph**  
research made easy

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## Contents

### Introduction

About this document.....	68
Typographic conventions.....	68
Copyright protection.....	68
Safety instructions in the official languages of the European Union .....	68

### General product information

Directives applied, product certification.....	69
California Residents .....	69
Customs declaration .....	69
Residual risk .....	69
Intended use.....	69
Compliant use.....	69
Reasonably foreseeable misuse .....	70
Transportation.....	70
Storage.....	70
Acclimatization .....	70
Permissible ambient conditions .....	70

### Safety

General safety instructions.....	71
Electrical safety.....	71
Data security .....	71
Operational safety .....	71
Occupational safety .....	72
Personal protective equipment (PPE).....	72
Environmental protection.....	72
Biohazard.....	72
Other regulations.....	73

### Device description

Mechanical design.....	74
Main Switch .....	75
Emergency off.....	75
Control panel .....	76
Interfaces .....	77
Sensor/peripherals connections .....	77

Cooling water supply.....	78
EASY LOCK flask clamping system .....	78
Mechanical design.....	79
Work positions .....	79
Handling .....	79
Clamping pressure adjustment .....	81
Flask support .....	81
Mechanical design .....	81
Handling .....	82
Heating Bath .....	83
Moving the heating bath manually .....	83
Automatic movement.....	83
Heating bath liquids.....	84
Automatic water refill heating bath .....	84
Safety features .....	85
Guard hood.....	85
Overpressure cut-out.....	86
Heating bath overheat protection .....	86
Heating bath filling level monitoring .....	86

## Commissioning

Positioning the device .....	87
Power supply .....	87
Switching the device on/off .....	88
Insert/remove the evaporation flask.....	88
Filling/emptying the heating bath.....	89
Bath liquid selector .....	89
Change the heating bath liquid .....	91

## Operation

Turn device control on/off.....	92
Structure of the user interface .....	92
Buttons and symbols .....	93
Numeric on-screen keyboard .....	94
Alphanumeric on-screen keyboard.....	94
[Settings] Menu .....	95
Security settings.....	96
Water timeout.....	97
Password / Screen lock.....	97
Temperature .....	98
Pressure.....	98

Bath liquid.....	98
Date/Time.....	99
Light inside Guard Hood .....	99
Glassware .....	99
Show Checklist .....	100
Show Checklist at Startup? .....	100
Language .....	100
Show System information.....	100
Datalog Settings .....	101
Show Datalog Series.....	101
Error history.....	101
Vacuum control mode .....	102
SETpressure mode.....	102
AUTOfaccurate mode.....	103
Distillation under atmospheric pressure .....	105
Set rotation speed .....	107
Recommended rotation speed.....	107
Set the heating bath temperature .....	108
Timer function .....	108
Favorites and Profiles.....	109
Favorites mode .....	109
Create/edit favorites .....	109
Load favorite.....	110
Live view.....	110
Create/edit process profiles.....	111
Load process profile .....	112
Data logging.....	113
Activating/deactivating data logging.....	114
Setting the logging rate .....	114
Logging capacity .....	114
Data list .....	114
Deleting measurement series.....	115
Saving measurement series on a USB flash drive .....	115
Data output .....	115
Troubleshooting.....	116

## Troubleshooting

Troubleshooting.....	116
----------------------	-----

## Appendix

Technical specifications .....	118
Scope of delivery .....	119
Accessories .....	119
Device service.....	120
General cleaning instructions.....	120
Cleaning and care instructions for heating baths .....	120
Repairs - Return of equipment.....	120
Maintenance.....	121
Disposal.....	121
Warranty statement.....	122
Contact details.....	122
Certificate of decontamination.....	123

### About this document

This operating instructions manual describes the features and operation of Hei-VAP Industrial large-scale evaporators. The operating instructions are an integral part of the described device!

### Typographic conventions

This document uses the following symbols, signal words, and highlights:

Symbol	Signal word / explanation
	Warning symbols in combination with a signal word indicate dangers: <b>DANGER</b>
	Indicates an immediate dangerous situation. Failure to respect the indications will result in death or serious injury.
	<b>WARNING</b>
	Indicates a potential danger. Failure to respect the indications will result in serious injuries.
	<b>CAUTION</b>
	Indicates a potential hazard which, if not avoided, damage to property and minor to moderate injuries can occur.
	Mandatory signs are used to indicate important and useful information on handling a product. This information is used to ensure operational safety and to maintain the value of the product.
→	The arrow indicates specific instructions to be followed to ensure operational safety when handling the product.

### Copyright protection

This document is protected by copyright and is intended for use by the purchaser of the product only.

No transfer to third parties, reproduction in any form, including excerpts, and by any means, as well as utilization and/or disclosure of the contents is permitted without the prior written consent of Heidolph Scientific Products GmbH. Any violation is subject to compensation for damage.

### Safety instructions in the official languages of the European Union

A summary of all safety instructions in the official languages of the European Union can be found in our Safety Guide for the product group Hei-VAP (shakers and mixers, Ref. 01-005-006-67). This document is available for download on our homepage in the most up-to-date version.

## Directives applied, product certification

### CE marking

The device meets all requirements of the following directives:

- European Machinery Directive, 2006/42/EC
- EMC Directive, 2004/108/EC

### NRTL Certification

The device has been tested in accordance with the following standards:

- UL 61010-1:2012/R:2016-04  
CAN/CSA-C22.2 No. 61010-1:2012/U2:2016-04
- UL 61010-2-010:2015  
CAN/CSA-C22.2 No. 61010-2-010:2015
- UL 61010-2-051:2015  
CAN/CSA-C22.2 No. 61010-2-051:2015

## California Residents

Important information for California residents regarding Prop 65. Please visit [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov) for more information.

## Customs declaration

Rotary evaporators, as combustion and cleaning equipment, may be subject to notification to a competent customs authority in the country of destination.

The evaluation of a customer obligation to notify and, if applicable, the notification to a competent customs authority in the country of destination is solely the responsibility of the user!

## Residual risk

The device has been designed and manufactured in accordance with the state-of-the-art standards at the time of development and the recognized safety regulations. During mounting and use, as well as during maintenance, repair and cleaning work, there are nevertheless certain residual risks associated with the described device.

These are identified and described at the appropriate points in this document.

## Intended use

The described device was designed by the manufacturer for the separation of substances, powder drying, concentration, crystallization of substances, and recycling of solvents under vacuum.

Due to its design, the device in its delivery condition may only be used in analytical processes or in laboratory-like conditions in the food, cosmetics, and pharmaceutical industries as well as other comparable industries that manufacture products intended for consumption by humans or animals, or for use on humans or animals.

Any other use of this device is not considered as intended!

## Compliant use

The user is generally responsible for evaluating the conformity of his application and, if necessary, for taking additional measures.

## **Reasonably foreseeable misuse**

Additional measures may be necessary for use under conditions or for purposes deviating from the intended use, and/or specific guidelines and safety regulations must be observed. Corresponding requirements must be evaluated and implemented by the operator in each individual case.

Compliance with and implementation of all relevant directives and safety measures for the respective field of application is the sole responsibility of the operator.

All risks resulting from improper use are borne solely by the operator.

The device may only be operated by authorized and instructed personnel. Training and qualification of the operating personnel as well as ensuring that the device is handled responsibly are the sole responsibility of the operator!

## **Transportation**

During transportation, avoid severe shocks and mechanical stresses that can cause damage to the device. Keep the original packaging in a dry and protected place for later use.

## **Storage**

Always store the device in its original packaging. To protect against damage and excessive material aging, store the device in an environment that is as dry, temperature-stable and dust-free as possible.

## **Acclimatization**

After each transportation and after storage under critical climatic conditions (e.g. high temperature difference between inside and outside), allow the device to acclimatize at room temperature for at least two hours to prevent possible damage due to condensation before commissioning it in the place of use. If necessary, extend the acclimatization phase if the temperature differences are very high.

Make all supply connections (power supply, tubing) only after the device has been acclimatized!

## **Permissible ambient conditions**

The device is designed for indoor use only. The device is **NOT** suitable for outdoor use! The device is **NOT** suitable for use in potentially explosive areas!

When used in corrosive atmospheres, the service life of the device may be reduced depending on the concentration, duration and frequency of exposure.

### General safety instructions

Before commissioning and using the device, familiarize yourself with all the safety regulations and occupational safety guidelines applicable at the place of use and observe them at all times.

Only operate the device if it is in faultless technical condition. In particular, ensure that there is no visible damage on the device itself and, where applicable, on connected devices or the supply connections.

If there is missing or misleading information on the device or regarding occupational safety, contact the responsible safety specialist or our technical service.

Only use the device in accordance with the regulations for intended use.

### Electrical safety

Before connecting the device to the power supply, ensure that the voltage indicated on the rating plate matches the specifications of the local power utility company.

Ensure that the power supply circuit provided is protected by means of a residual-current device (RCD).

Always use the power supply cord provided with the device. All installation work for the power supply of the device may only be carried out by an authorized electrician or by the technical service of Heidolph Scientific Products GmbH!

Prior to use, check that the device and the power supply cord are free of visible damage.

Have repairs and/or maintenance work on the device carried out exclusively by an authorized electrician or by the technical service department of Heidolph Scientific Products GmbH.

Always switch off and disconnect the device from the power supply, preventing reconnection, before carrying out maintenance work, cleaning, or repairs.

### Data security

The user is responsible for ensuring data security when transferring data between the described device and other devices.

Only use secure networks for the data transfer and avoid use of critical infrastructure.

Only use high-quality shielded data cables for the data transfer.

For data transfer via a USB B connection, an industrial standard USB hub should be preferably used to ensure the most stable connection possible.

### Operational safety

Operate the device under a closed ventilated fume hood when working with potentially hazardous substances (see EN 14175 and DIN 12924).

Do not make any unauthorized changes or modifications to the device!

Only use genuine spare parts and accessories, or those expressly approved by the manufacturer!

Rectify malfunctions or faults on the device immediately.

Switch off and disconnect the device from the power supply, preventing reconnection, if it is not possible to eliminate the malfunction or rectify the fault immediately.

Observe all relevant general and safety instructions for the connected peripheral devices (observe the supplied documentation!).

Observe all other applicable regulations such as laboratory and workplace guidelines, recognized safety technology rules and special local regulations.

## Occupational safety

Always use the prescribed personal protective equipment (PPE) such as protective clothing, safety goggles, protective gloves, safety shoes, etc.

Do not operate any other devices in the immediate vicinity of the device ...

- which can generate electromagnetic fields in the frequency range between  $9 \times 10^3$  Hz to  $3 \times 10^{11}$  Hz,
- which generate emission or radiation sources in the frequency range  $3 \times 10^{11}$  Hz to  $3 \times 10^{15}$  Hz (in the optical spectral range wavelengths from 1,000 µm to 0,1 µm),
- which generate ultrasonic or ionizing waves.

Do not operate the unit when adiabatic compression or shock waves may occur (pressure wave ignition).

Do not use substances that could release energy in an uncontrolled way and cause a pressure increase (exothermic reaction, spontaneous ignition of dusts).

Do not process hard, brittle materials such as stones, soil samples, etc., that could destroy the evaporation flask.

Only use heating bath media that guarantee sufficient heat transfer.

Do not operate the device with overpressure.

Do not expose the glass components to a pressure difference of more than 2 bar.

Ensure that the coolant overpressure does not exceed a level of 2 bar.

Ensure that the flow velocity does not exceed 1 m/s when aspirating liquids with flammable components (electrostatic charge, danger of ignition!).

Avoid the formation of explosion group IIC gases and potentially explosive distillation residues.

## Personal protective equipment (PPE)

The operating company must determine and provide the necessary PPE, depending on the respective application and the media and chemicals used.

The required instruction of the personnel is solely within the operating company's responsibility.

## Environmental protection

When processing environmentally hazardous substances, take appropriate measures to avoid hazards to the environment.

The evaluation of appropriate measures such as the marking of a hazardous area, their implementation, and the training of the relevant personnel is the sole responsibility of the operator!

## Biohazard

When processing biohazardous substances, take appropriate measures to prevent hazards to persons and the environment, including:

- Instruction of the personnel regarding the necessary safety measures.
- Provision of personal protective equipment (PPE) and instruction of the personnel in its use.
- Marking of the device with the biohazard warning symbol.

The evaluation of appropriate measures such as the marking of a hazardous area, their implementation, and the training of the relevant personnel is the sole responsibility of the operator!

## Other regulations

In addition to the notes and instructions in this document, observe all other applicable regulations such as laboratory and workplace guidelines, hazardous substances ordinances, recognized rules of safety engineering and occupational medicine as well as particular local regulations!

Noncompliance will invalidate any warranty claims against Heidolph Scientific Products GmbH.

The operator is solely liable for all damage resulting from unauthorized changes or modifications to the device, from the use of unapproved or non-genuine spare parts and accessories, or from disregarding the safety instructions and hazard warnings or the manufacturer's instructions!

## Mechanical design

The following figure shows an example configuration with glassware RC and optionally available base cart.

### Index

<b>1</b>	Base cart (option)
<b>2</b>	Control panel
<b>3</b>	Evaporation flask
<b>4</b>	Guard hood
<b>5</b>	Emergency off switch
<b>6</b>	Rotation axis locking
<b>7</b>	Glassware RC, including condenser and receiving flask
<b>8</b>	Condenser protection with PMMA disc

### Main Switch

The main switch for switching the device on and off is located on the rear of the device.

Turn the main switch clockwise to the ON position to switch on the device.

Note that the mains voltage is also present at the device when it is switched off (main switch in the Off position). In order to de-energize the device for maintenance or repair work, the device must be disconnected from the mains!

### Emergency off

The emergency off button of the device is located on the front, above the control panel.

By pressing the emergency off button, all device functions are switched off and the running process is stopped, but the mains voltage is still present! The device switches to the safe state (rotation and heating off, heating bath moves to the lower end position), but there is still mains voltage!

In the event of an emergency off, a warning signal sounds for approximately 5 seconds.

The emergency off button engages when pressed. Correct the cause of the error and proceed as follows to unlock the emergency off switch:

- Turn off the device
- Unlock the emergency off button, turning it slightly to the left.
- Switch the device back on after a wait of at least 10 seconds

As soon as the cause of the error has been corrected and the emergency off button has been released, the device is ready for operation again.

Before switching the device back on after an emergency off, it is essential to correct the cause(s) of the fault or disturbance.

Use the emergency off button only to switch off the device in emergency situations or in the event of a fault.

Always use the main switch to switch off the device during normal operation!

## Control panel

### CAUTION

Pointed or sharp objects can damage the surface of the display on the control panel!

To operate the device, touch the display only with your fingertips or a special stylus with a soft-touch tip.

The control elements of the central display are touch-sensitive (touch panel) and can be operated by finger or with a suitable stylus.

### Index

- |           |   |
|-----------|---|
| <b>1</b>  | Menu button   |
| <b>2</b>  | Mode selection buttons  |
| <b>3</b>  | [Start/Stop] button   |
| <b>4</b>  | [Start/Stop Timer] button   |
| <b>5</b>  | Heating bath lift up/down buttons   |
| <b>6</b>  | [Process Parameters] area with mode button, set/actual display and start/stop button for individual processes |
| <b>7</b>  | [General Process Stop] button: stop all processes simultaneously  |
| <b>8</b>  | [Warm-up] button  |
| <b>9</b>  | [General process start] button: start all processes simultaneously  |
| <b>10</b> | [Vacuum override] control   |
| <b>11</b> | On/off switch   |

## Interfaces

The device interfaces are located on the bottom of the control panel:

- 1 x USB 2.0, max. 500 mA
- 1x LAN RJ-45 (without function in the current device version)

### **WARNING**

If the interfaces are assigned too high voltages and are not sufficiently insulated, metallic parts such as the housing can be live in the event of a fault.

Safely isolate ultra-low-voltage inputs and outputs via 25 V AC or 60 V DC according to EN 61140, or by double or reinforced insulation according to EN 60730-1 or DIN 60950-1.

Use only shielded connection cables. Connect the shield to the connector housing.


### **CAUTION**

Make sure that the device interfaces on the bottom of the display are always properly covered!

Penetrating liquid can cause interface malfunctions and short-circuiting.

## Sensor/peripherals connections

The connections for the sensors and the peripheral devices are located on the rear panel of the device (further information, see mounting instructions O1-O01-O09-12):

Connections	
1	Vacuum sensor
2	Vacuum valve
3	Venting valve
4	Heating bath water valve
5	Vacuum pump
6	Domestic water inlet automatic water refill
7	In: cooling liquid from recirculation chiller
8	Out: cooling liquid to recirculation chiller
9	Emergency off sensor Distimatic Pro (optional)
10	Temperature sensor PT (cooling liquid recirculation chiller)

## Cooling water supply


### CAUTION

Keep the supply and return connections closed if no external condenser is connected!

In order to not exceed the maximum permissible pressure in the chiller (2 bar) during operation, open the supply and return connections when an external condenser is connected!

The connections for the cooling water supply are located on the left-hand side panel of the device, above the receiving flask(s). The feed of the condenser is made via the connections [2] and [4] (main cooling, cooling of the medium). Connections [1] and [3] are equipped with shut-off valves and can be used for optional units:

#### Connections

- |          |   |
|----------|---|
| <b>1</b> | Cooling liquid supply optional units                |
| <b>2</b> | Cooling liquid return from condenser (main cooling) |
| <b>3</b> | Cooling liquid return optional units                |
| <b>4</b> | Cooling liquid supply condenser (main cooling)      |

## EASY LOCK flask clamping system

The device is equipped with the patented EASY LOCK flask clamping system.

### CAUTION

Incorrect or careless handling can damage the flask clamping system and/or the flask inserted!

Before inserting a flask for the first time, familiarize yourself with the instructions for handling the EASY LOCK flask clamping system.

Before inserting a flask for the first time, familiarize yourself with the instructions for handling the flask support (see "Flask support" on page 81).

Unsuitable flasks can damage the flask clamping system and/or the flask inserted!

The EASY LOCK flask clamping system is only suitable for flasks with an opening diameter of 150 mm. The use of flasks with a smaller opening diameter (e.g. LR 20) is not permitted, as these cannot be fixed in a positive manner due to the geometric properties.

## Mechanical design

Components	
1	Clamping lever
2	Clamping lever lock
3	Outer ring
4	Sliding segments for fixing the flask

## Work positions

The following figure shows the different working positions of the flask clamping system:

- "closed-locked" position [1]: In this position, the evaporation flask is clamped in the EASY LOCK, the evaporator can be operated normally.
- "unlocked" position [2]: The internal sliding segments of the EASY LOCK are open enough that the evaporation flask inserted can be rotated and positioned for a change.
- "open" position [3]: The evaporation flask can be inserted or removed. Furthermore, the sliding segments and the outer ring can be removed.

## Handling

Before opening the EASY LOCK, the rotation axis of the device must be locked using the latching rotary knob!

- Pull out the latching rotary knob and turn it up to the lock position to lock the rotation axis.
- Depending on the rotation position, it may not be possible to engage the locking knob directly in the locking position! In this case, turn the outer ring of the EASY LOCK and the latching rotary knob at the same time until it can be turned to the locking position and locks into place.
- Turn the latching rotary knob into the unlock position to unlock the rotation axis.


After locking the rotational axis, the following message appears on the device display:


A process start is only possible when the rotational axis is unlocked and the warning message on the display is confirmed with [OK]!

- To open the EASY LOCK, first press in the clamping lever lock to unlock the clamping lever.
- Pull the unlocked clamping lever up to the stop until it is held in the open position (2) by the pin on the threaded bolt of the pretensioning mechanism.  
The sliding segments are released, the flask can be moved but not removed!
- Grip the clamping lever with your hand and turn the outer ring of the EASY LOCK downwards up to the stop (3). When doing so, the clamping lever must not be pressed against the outer ring.  
The EASY LOCK is now fully open and the flask and the PTFE seal can be inserted or removed.

- Insert a suitable evaporation flask.
- Grip the clamping lever with your hand and use it to pull the outer ring of the EASY LOCK upwards up to the stop.  
The sliding segments are pre-tensioned, the flask can be moved, but is securely held.
- Then press the clamping lever against the outer ring with normal manual force until the clamping lever lock audibly engages.

### Clamping pressure adjustment


Due to industrial manufacturing tolerances, the contact pressure of the EASY LOCK must be adjusted individually for each evaporation flask.

To do this, proceed as follows:

- Lock the rotational axis.
- Unlock and open the EASY LOCK clamping lever to the stop position.
- Insert the set screw completely, turning clockwise (marked in the following figures): In this position, minimum contact pressure is applied to the flask.
- Insert the evaporation flask and lock the EASY LOCK.
- Check if the clamped flask can be moved in the locked EASY LOCK.

#### If the flask still rotates:

- Unlock and open the EASY LOCK clamping lever to the stop position.
- Pull back the inserted flask.
- Turn the set screw one to two turns counter-clockwise to increase the contact pressure step by step.
- Insert the evaporation flask and lock the EASY LOCK.
- Check if the clamped flask can be moved in the locked EASY LOCK.
- Repeat the procedure until the clamping lever closes and locks against moderate resistance and the flask can no longer be moved by hand in the closed EASY LOCK

#### CAUTION

Excessive clamping/closing pressure can damage the EASY LOCK and/or the flask inserted!

Always follow the instructions for the correct adjustment of the contact pressure in this section!

### Flask support

The flask support facilitates the insertion of the evaporation flask and prevents the flask from falling into the heating bath during removal.

### Mechanical design

The flask support consists of the following mechanical components (see also the following figure):

- Two guide bars, fixed to the chassis
- Two support bars with support element, movable on the two guide bars
- Supported locking lever with handle and defined locking positions

### Components

<b>1</b>	Guide bars
<b>2</b>	Support bars with flask support element
<b>3</b>	Locking lever

## Handling

The two support bars of the flask support are connected with the locking lever so that both can be moved parallel along the guide bars to and from the EASY LOCK.

The distance between the support bars is defined by the locking positions for three flask sizes.

To set the appropriate distance, do the following:

- Uncouple the support bars by pulling up on the locking lever on the handle.
- Specify the appropriate locking position (refer to the volume mark on the evaporation flask!).
- Move the two support bars so that the locking lever engages in the appropriate locking position when moving it down.
- Using the locking lever, re-couple the support bars.
- Ensure that the support bars are securely coupled. To do this, move the unit along the guide bars.

The evaporation flask is held by the support bars in such a way that it can be rotated in all axles.

Due to the mechanical design of the flask supports, the evaporation flask cannot fall into the heating bath or roll from the flask support when the locking lever is correctly set and engaged.

For normal operation, the two support bars must be decoupled and moved outwards to their respective end position. The end position of the support bars is monitored electronically.

As long as the support bars are not in their respective end position, the following message appears on the display:

A process start is only possible when the support bars are in their respective end position and the warning message on the display has been confirmed with [OK]!

### Heating Bath

The heating bath is moved up or down with a motor lift. The control is carried out via the device software and the control panel of the device.

In normal operation, the heating bath is moved to the upper end position in order to dip the flask into the heating bath liquid.

To replace the evaporation flask, to fill the heating bath and to carry out cleaning, maintenance or repair work, etc., the heating bath is moved to the lower end position.

The heating bath lift only works if the level sensor is correctly mounted in the heating bath.

### Moving the heating bath manually

The [lift DOWN/UP] arrow keys can be used to manually move the heating bath lift. The buttons are available on the start window after the device is switched on and initialized:

- Touch the [lift DOWN] button to move the heating bath down manually.
- Touch the [lift UP] button to move the heating bath up manually.

The movement stops as soon as you release the button or automatically when the upper or lower end position is reached.

#### CAUTION

When the evaporation flask is immersed, the already filled heating bath liquid is displaced and the liquid level in the heating bath rises!

Especially when using oil, the device and the surrounding floor surface can be heavily contaminated by overflowing heating bath liquid.

When filling the heating bath, make sure that the maximum filling quantity is not exceeded.

### Automatic movement

The heating bath is automatically moved up or down when one of the following device functions is performed:

- Warm Up: when the heating process is started, the heating bath is automatically moved to the lower end position - the flask is lifted out of the heating bath.
- Start All: when starting all processes in automatic mode, the heating bath is automatically moved to the upper end position - the flask is immersed in the heating bath liquid.
- Stop All: when stopping all processes, the heating bath is automatically moved to the lower end position - the flask is lifted out of the heating bath.
- Emergency off: when the emergency off button is pressed, the heating bath is automatically moved to the lower end position - the flask is lifted out of the heating bath.
- Timer:Lift: after a set time period in timer mode, the heating bath is automatically moved to the lower end position - the flask is lifted out of the heating bath.

## Heating bath liquids

### DANGER

When using oils with a flash point < 285 °C as a heating bath liquid, uncontrolled thermal reactions can occur at high heating bath temperatures and an explosive atmosphere can arise: Only use oils with a flash point > 285 °C as heating bath liquid.

From a heating bath temperature of approx. 75 °C – 80 °C, the use of polyethylene glycol or oil is strongly recommended.

Reference values heating bath liquid:

- Tap water; temperature range 20 °C – 100 °C
- Polyethylene glycol; temperature range 20 °C – 180 °C
- Silicone oil, low viscosity (max. 40 cP), flash point > 285 °C; temperature range 20 °C – 170 °C

## Automatic water refill heating bath

This function is activated when the heating bath liquid selector is set to [Water] and the corresponding setting has been made: see "Bath liquid selector" on page 89.

### Functional principle

Due to the heat development in the heating bath, the fluid level gradually decreases due to evaporation.

In order to ensure a constant immersion depth of the flask and thus a continuous evaporation process, the heating bath is automatically lifted up step by step toward the upper end position. This lifting movement is controlled by the level sensor.

### Operating behavior

As soon as the heating bath reaches the upper end position (minimum level), water is automatically refilled until the maximum level in the heating bath is reached again.

### Monitoring

When refilling the heating bath via the automatic water refill, the filling time is limited to 30 seconds for safety reasons. To deactivate the function or to personalize the filling time, open the [Settings] menu and edit [Water Timeout] (see section "Water timeout" on page 97).

### CAUTION

Possible overflow when refilling the heating bath! The monitoring of the refill duration during the automatic water refill should only be deactivated in exceptional cases.

If the maximum level during an active automatic water refill is not reached within the defined time period, the automatic refill is shutdown and the device switches to the safe operating state.

The following message appears on the display:

### Error

Water Timeout exceeded.  
Autom. heating bath filling stopped!  
Check heating bath level sensor and  
water inlet.  
Restart unit via the main switch on  
the backside.

- Restart the device. Check the water supply from the house connection to the supply above the heating bath for correct connection and continuity as well as the level sensor for correct function.

If the maximum level is reached within the defined time period, the automatic refill is switched off normally and the device is ready for use.

## Safety features

### Guard hood

The device is equipped with a guard hood that completely covers the evaporation flask and the heating bath, see section "Mechanical design" on page 74.

---

#### **WARNING**

Even when the guard hood is closed, there is a risk that the heating bath liquid will spill out of the heating bath at high rotational speeds. There is a risk of scalding and/or contamination.

Also refer to the recommendations in section "Set rotation speed" on page 107.

#### **CAUTION**

Always use the handle provided to open and close the guard hood.

When closing the guard hood, keep sufficient distance to the closing edges with your free hand.

---

The guard hood is equipped with an interior lighting which allows a better view of the ongoing distillation process. The guard hood fulfills the following safety functions:

- Basic protection against touching the rotary/heated evaporation flask
- Basic protection against direct contact with the (heated) bath liquid
- Basic protection against flying glass parts and spritating liquid in the case of implosion or thermal reaction

An error message appears on the display when the guard hood is opened during operation. At the same time, the flask rotation is stopped and all running processes are stopped. As soon as the guard hood is fully closed again, the flask rotation starts again and the stopped processes continue.

Due to design and function, the guard hood is not designed to tightly seal the heating bath area!

### Overpressure cut-out

The device is equipped with an electronic overpressure monitoring system. As soon as the integrated vacuum sensor detects an impermissible pressure level, an error message appears on the display; at the same time, the device is switched off and put into an operationally safe state.

The vacuum limit is set to 1,200 mbar at the factory. Depending on the prevailing ambient conditions, this limit value can be adjusted in a range of 900 – 1,400 mbar.

The overpressure limit is set in the [Settings] – Safety Settings menu, see "Security settings" on page 96.

#### PRACTICAL TIP

The higher the application location is above sea level, the lower the vacuum limit value:

- Measure the ambient pressure on the ventilated system and at high air pressure.
- Add about 50 mbar to the reading.

### Heating bath overheat protection

The device is equipped with an autonomous overheat protection function that prevents the heating bath from overheating. When reaching the maximum heating bath temperature, the heating function is switched off by a mechanical switch

The overheat protection must be reset manually:

- Make sure that the heating bath is in the lowest position and allow the heating bath liquid to cool down.
- Remove the safety assembly cover on the right side of the rotary evaporator housing, see figure.
- Operate the reset button with a suitable tool (blunt tip, insulated or non-conductive material).
- Mount the safety assembly cover on the right side of the rotary evaporator housing.

### Heating bath filling level monitoring

The heating bath filling level is monitored by a permanently installed level sensor (float switch). The movement of the heating bath lift is stopped when the level sensor is triggered. Level sensor limit values:

- minimum liquid level in the heating bath: 50 mm from the bottom of the heating bath
- maximum liquid level in the heating bath: 145 mm from the bottom of the heating bath

#### Operating behavior with automatic refilling function

For safety reasons, the filling time is limited to five (5) minutes when initially filling or refilling the heating bath via the automatic water refill function. This value cannot be changed by the user!

If the maximum filling level is not reached within the defined period of time, the automatic filling function is turned off and the device switches to a safe operating state.

- In this case, check the water supply from the house connection to the supply above the heating bath for correct connection and continuity as well as the level sensor for correct function.

If the maximum level is reached within the defined time period, the automatic refill is switched off normally and the device is ready for use.

## Positioning the device

Ensuring proper system and/or device mounting, including all accessory parts, is the sole responsibility of the operating company.

If required, Heidolph offers a professional system and/or device mounting service, including all supplied components and including commissioning. The associated mounting instructions (Ref. O1-001-009-12) describe all steps for proper system and/or device mounting in detail!

---

### CAUTION

Due to improper mounting and/or installation of the system/device, as well as unauthorized changes to the mounting, there is a risk of direct and indirect damage to property!

If necessary (system mounting and installation carried out by the operator or third parties, necessary changes to the existing mounting), contact the manufacturer's technical service (see section "Contact details" on page 122).

Place the device on a clean, stable, level and horizontal surface.

Before switching on, make sure that the device is sufficiently stable and use all existing locking devices on the base cart.

---

The operator is solely liable for all damage resulting from unauthorized changes or modifications to the device, from the use of unapproved or non-genuine spare parts and accessories, or from disregarding the safety instructions and hazard warnings or the manufacturer's instructions! Any warranty claim against Heidolph Scientific Products GmbH shall lapse in such a case

## Power supply

---

### DANGER

Observe the instructions in section "General safety instructions" on page 71, in particular Electrical safety.

---

## Switching the device on/off

### CAUTION

Before each operation of the device, check the caps and connecting elements of all glass components for correct and tight fit!

At loose/leaky connections, ambient air is sucked in during the vacuum build-up and the required vacuum pressure cannot be reached in a stable manner!

In the event of excess pressure, there is a risk that liquid will escape from the system! The device and surrounding areas can be severely contaminated by leaking media!

Tighten loose caps/fasteners by hand or correct the seating of the retaining clips.

To switch on and off, use the main switch on the rear of the device, see section "Main Switch" on page 75.

## Insert/remove the evaporation flask

### CAUTION

In the event of improper handling, the evaporation flask that is to be used may fall to the ground.

- Glass breakage!
- Risk of injury/contamination due to leaking media!

Use only the evaporation flasks approved by the manufacturer.

Pay particular attention to all information on handling the "EASY LOCK flask clamping system" on page 78.

### WARNING

The smallest damage to the evaporation flask can lead to glass breakage when pressurized.

- Glass breakage!
- Risk of injury/contamination due to leaking media!

Before each operation of the device, check the evaporation flask for visible damage.

## Insert the evaporative flask

To insert the evaporation flask, follow these steps:

- Stop all running processes (if applicable).
- Move the heating bath to the lower end position ("Heating Bath" on page 83).
- Stop the rotation movement of the flask (if applicable).
- Open the guard hood.
- Lock the rotation axis and fully open the flask clamping system.
- Position the flask support over the heating bath according to the flask size.
- Align the evaporation flask on the flask support so that the opening can be inserted into the open EASY LOCK.
- For new evaporation flasks, adjust the contact pressure of the EASY LOCK.
- Close the EASY LOCK.
- Make sure that the EASY LOCK clamping lever is properly engaged.

### Removing the evaporation flask

To remove the evaporation flask, follow these steps:

- Stop all running processes (if applicable).
- Move the heating bath to the lower end position ("Heating Bath" on page 83).
- Stop the rotation movement of the flask (if applicable).
- Lock the rotational axis.
- Open the guard hood.

---

#### CAUTION

The heating bath is designed for a maximum operating temperature of 210°C (using silicone oil). From a liquid temperature of 50 °C, there is a risk of injury in the event of contact.

Allow heated liquid to cool below 50 °C before removing the flask.

Use appropriate personal protective equipment (heat-resistant gloves, eye protection, safety clothing) to remove the flask.

---

- Position the flask support according to the flask size.
- Fully open the flask clamping system.
- Using the flask support, pull the evaporation flask out of the open EASY LOCK and rotate it so that the opening is vertical up.
- Removing the evaporation flask.

### Filling/emptying the heating bath

Before the initial filling or refilling of the heating bath, the preselection water or oil has to be made via the [Bath liquid] selector.

#### Bath liquid selector

Before filling the heating bath for the first time and when changing the bath liquid, the bath liquid selector on the back of the device must be moved to the required position.

The bath liquid selector setting affects the following device functions:

- Process control (automatic mode)
- Heating bath level sensor
- Automatic water refill heating bath

The selection can only be changed if the device is switched off and the heating bath liquid has cooled down to at least 50 °C!

- Slide the bath liquid selector to the oil or water position.
- [Water] position
  - Maximum temperature setting: 90 °C
  - Automatic water refill active
  - Heating bath level sensor active
- [Oil] position
  - Maximum temperature setting: 180 °C
  - Automatic water refill inactive
  - Heating bath level sensor deactivated
- The setting of the bath liquid selector must be confirmed in the [Settings] menu, see "Bath liquid" on page 98.

### Fill the heating bath

The heating bath can be filled manually in any case and automatically when tap water is used. Before automatic filling, make sure that the device is correctly connected to the house connection. Observe the instructions in section "Operating behavior with automatic refilling function" on page 86.

---

#### CAUTION

If the filling level is too high, there is a risk that the heating bath liquid will overflow when the flask is immersed, possibly contaminating severely the surrounding area!

Observe the mark for maximum filling level in the heating bath when manually filling!

- Set the bath liquid selector on the rear of the device to the correct position, see section "Bath liquid selector" on page 89.
- Close the drain valve on the back of the device (handle in 90 ° – position to the drain valve).
- Ensure that the evaporation flask is correctly mounted.
- Move the heating bath to the upper end position.
- Fill the heating bath with the flask immersed, until the liquid level reaches a height of 40 mm below the top edge (maximum filling height) or open the valves from the house connection for filling the heating bath via the filling system of the device.

### Empty the heating bath

---

#### CAUTION

The heating bath is designed for a maximum operating temperature of 210°C (using silicone oil). From a liquid temperature of 50 °C, there is a risk of injury in the event of contact.

Allow heated liquid to cool below 50 °C before removing the flask.

Use appropriate personal protective equipment (heat-resistant gloves, eye protection, safety clothing) to remove the flask.

- Ensure that a suitable drain tube is attached to the drain valve.
- Move the heating bath to the lower end position ("Moving the heating bath manually" on page 83).
- Open the drain valve on the back of the device (handle parallel to the drain valve).
- Clean the heating bath and the drain tube between the heating bath and the drain cock.
- Close the drain valve on the back of the device!

### Change the heating bath liquid

To change the heating bath liquid, proceed as described in the preceding section "Removing the evaporation flask" on page 89:

- Emptying the heating bath, see "Empty the heating bath" on page 90.
- Filling the heating bath, see "Fill the heating bath" on page 90.
- Observe the information given in section "Bath liquid selector" on page 89!

When changing the bath liquid, the heating bath and the drain tube between the heating bath and the drain valve must be thoroughly cleaned and dried.

In particular, when changing from water to oil, make sure that the heating bath temperature is only slowly and gradually augmented during restart (recommendation: 10 °C), so that remaining water in the system can evaporate until the process temperature is reached.


### Turn device control on/off

Make sure that the device is switched on (main switch ON).

To enable/disable the device control, touch the on/off button in the lower right corner of the control panel (see "Control panel" on page 76).

First, the shown dialog box appears on the device display.

- Confirm the prompt with [Yes], to deactivate the device control.
- Confirm the prompt with [No], to activate the device control.

### Structure of the user interface

After the device control has been initialized, the user interface start page appears on the display.

The following options are available:

- Menu [Settings]: To open the menu, touch the [Menu] button in the upper left corner.
- Selection of the operating mode
  - Touch the [Set Pressure] button to operate the device at a specific vacuum pressure, see "SETpressure mode" on page 102.
  - Touch the [Auto Accurate] button to activate the AUTOaccurate function and operate the system with automatic vacuum control, see "AUTOaccurate mode" on page 103.

The software keys/buttons of the active function are highlighted in orange in all menu windows.

## Buttons and symbols

The following table provides an overview of all buttons and icons on the control panel and in the menu navigation:

Symbol	Function
	ON/OFF button on the control panel
	[Menu] button
	Start/stop function
	Heating bath lift up / down
	Icon [function disabled]
	Confirm (OK), activate function / Icon [function activated]
	Open info dialog (checklist)
	Jump to the end (checklist)
	Arrow keys step-by-step back/forward
	[Delete] button
	Toggle keys numeric/alphanumeric characters
	Shift button upper case
	Insert space
	Add new favorite Add new row (Profile)
	Remove individual entry
	Remove entire list or operation
	Zoom in/out
	Save Favorite
	Save and apply Favorite

## Numeric on-screen keyboard

The numeric on-screen keyboard appears automatically when you touch a numeric entry field within a mask.

- Use the keyboard to enter a desired value within the range of values allowed in each case.
- Touch the [Cancel] button to return to the parent dialog without making any changes.
- Touch the [Cancel] button to cancel the last entry.
- Confirm your entry with [OK] to accept the new value. The on-screen keyboard closes.

1	2	3
4	5	6
7	8	9
.	0	
Cancel	☒	Ok

## Alphanumeric on-screen keyboard


The alphanumeric on-screen keyboard appears automatically when you touch an input field for alphanumeric values within a mask.

- Use the keyboard to enter a desired designation.
- Touch the [Cancel] button to return to the parent dialog without making any changes.
- Touch the [123] button to open the highlighted keyboard with numeric values and special characters.
- Touch the [Cancel] button to cancel the last entry.
- Confirm your entry with [OK] to accept the new value. The on-screen keyboard closes.

a	b	c	d	e	f
g	h	i	j	k	l
m	n	o	p	q	r
s	t	u	v	w	x
↑	,	.	-	y	z
Cancel	123	←	☒	Ok	

## [Settings] Menu

From the start window, touch the [Menu] button to open the main menu. Touch the [Settings] button to open the [Settings] menu. Here you can set various default settings for the operation of the device:

The menu contains several screens. Touch the display and swipe up or down to scroll through the menu. The operating functions are identical in all screens:

- Touch the [right] arrow key of an entry to open its submenu.
- Touch an entry field to open the on-screen keyboard and enter a desired value.
- Touch the [left] arrow key at the bottom of the screen to return to the [Settings] menu without any changes.
- Touch the [OK] button at the bottom right of the screen to confirm an entry.

Valid changes to the parameters will be accepted without any further security prompt! Inputs outside the valid value range are not accepted.

## Security settings

### CAUTION

Due to incorrect/insufficient safety settings there is a risk that processes will not be shut down or will be shut down with a time delay when an unacceptable pressure level is reached.

Only qualified personnel may adjust the safety settings of the device.

When required, please contact our technical service, see "Kontaktdaten" on page 61 section.

Safety Settings		
Overpressure	1400	mbar
Leakage Pressure	900	mbar <input checked="" type="checkbox"/>
Overheat protection	10	°C <input checked="" type="checkbox"/>

In the [Settings] menu, the conditions under which a running process is automatically stopped can be defined.

### Overpressure

- Touch the [Overpressure] data entry field and use the on-screen keyboard to enter the desired value.
  - Default [Overpressure]: 1200 mbar
  - Setting range [Overpressure]: 900 – 1400 mbar
- Confirm with [OK].

### Leakage pressure

The [Leakage Pressure] stop criterion must be manually activated/deactivated. Default setting: Inactive:

- Touch the [function disabled] icon of the entry. The icon changes to [function activated], the stop criterion is active and can be adjusted.
- Touch the [Leakage Pressure] data entry field and use the on-screen keyboard to enter the desired value.
  - Default [Leakage Pressure]: 900 mbar
  - Setting range [Leakage Pressure]: 1 – 1399 mbar
- Confirm with [OK].

### Overheat protection


The [Overheat Protection] stop criterion must be manually activated/deactivated. Default setting: Inactive:

- Touch the [Function disabled] icon of the entry. The icon changes to [function activated], the stop criterion is active and can be adjusted.
- Touch the [Overheat Protection] data entry field and use the on-screen keyboard to enter the desired value.
  - Default setting: 15 °C
  - Adjustment range: 10 – 25 °C in increments of 1 °C.
- Confirm with [OK].

## Water timeout

In this menu, the refill time for the [Automatic water refill] function can be adjusted and deactivated.

- Adjust refill time, [function activated] icon: Touch the entry field for seconds [s]:
- Enter the desired waiting time in a range of [1 s – 9999 s] using the on-screen keyboard.
- Deactivate function: Touch the [function activated] icon. The icon changes to [function disabled]:


## Password / Screen lock

The following options are available in this menu:

### Password

After a specified period of inactivity (see the following parameter [Lock screen after]), the device's panel is automatically locked. Set a password here to secure the device against unauthorized access.

- Touch the [function disabled] icon to open the password entry screen.
- Set a freely selectable password (at least eight characters) for sharing the lock screen: Use the on-screen keyboard to enter the desired password on the first line. Enter the password again in the second line and confirm with [OK].

If password protection is active, the vacuum controller is also deactivated in the locked state, regardless of the setting [Manual vacuum override active while screen locked], see below.

When password protection is active, the lock symbol is shown in orange on the display.

Incorrect password entries, regardless of frequency, have no effect on an ongoing process. If the password is lost, the device can be restarted and a new password set if necessary. However, when the device is restarted, the current process is aborted and cannot be continued after the device is restarted.

To disable password protection, touch the [function enabled] icon. The symbol changes to the [function deactivated] symbol, the password protection is inactive.

### Lock screen after

You can use this menu item to specify the time after which the display is automatically locked from the last operation.

- Touch the [Minutes] entry field.
- Use the on-screen keyboard to enter the desired wait time in a range of at least one (1) to a maximum of five (5) minutes.

To unlock the screen, drag the slider to the right.

### Manual vacuum override active while screen locked

As soon as the screen lock is active, the system stops responding to the vacuum control knob of the control panel in its normal state. This menu item allows you to deactivate the vacuum control knob lock.

- Touch the [function disabled] icon of the entry. The icon changes to [function activated]. In this state, the vacuum pressure can be adjusted using the vacuum control knob despite the screen lock being active.

Proceed in the same way to activate the vacuum control knob lock again. Touch the [function enabled] icon, icon changes to [function disabled].

### Temperature

You can use this menu item to select the unit for temperature displays.

- Touch the button of the desired unit (°C, °F, or K).

The active selection is highlighted in orange. All temperature values are displayed in the selected unit.

### Pressure

You can use this menu item to select the unit for pressure displays.


- Touch the button of the desired unit (mbar, hPa, or Torr).

The active selection is highlighted in orange. All pressure values are displayed in the selected unit.

### Bath liquid

The current setting for the bath liquid is displayed here.


- First check the position of the bath liquid selector on the back of the device (see section "Bath liquid selector" on page 89) and adjust it if necessary.
- Touch the [Change] button to switch between the [Water] and [Oil] setting.
- Note and confirm the shown advice with [OK]:

## Date/Time

This menu item allows you to adjust the format for displaying the system time.

- Touch the [right] arrow key in the entry. The [Date/ Time] input screen appears on the display.
- Set the date and time format by touching the assigned input field (active format highlighted in orange).
- Touch the [Date] and [Time] input fields one after the other and enter the desired values using the on-screen keyboard.

## Light inside Guard Hood

This menu item allows you to activate or deactivate the automatic activation of the interior lighting of the guard hood.

- Touch the [function disabled] icon of the entry. The icon changes to [function activated]. In this state, the interior lighting of the guard hood is automatically switched on when the device is switched on.


Follow the same procedure to switch off the automatic activation of the interior lighting again: Touch the [function enabled] icon, icon changes to [function disabled].

## Glassware

You can select the mounted glassware via this menu item.

- Touch the [right] arrow key in the entry. The [Glassware Assembly] selection screen appears on the display.
- Touch the entry of the mounted glassware (A, A2 or R) to load the corresponding system settings. The settings are immediately applied in the main menu and in the checklist.

### Show Checklist

You can use this menu item to call up a system checklist. You can use this checklist to check all the relevant components on the device structure: Basic safety instructions, installation of peripheral devices, glassware, sensor system and tubing, as well as process aid.


- Touch the arrow key [right] in the entry and confirm the system query [Show Checklist] in the following window with [OK] to load the checklist stored for the selected glassware.

### Show Checklist at Startup?

You can use this menu item to specify whether the checklist is loaded automatically when the device control is switched on (factory setting!), or manually as described above.


- Touch the [function enabled] icon of the entry. The icon changes to [function disabled]. In this state, the checklist is not loaded when the device control is switched on.
- Follow the same procedure to enable automatic loading of the checklist again. Touch the [function enabled] icon, icon changes to [function disabled].

### Language

You can use this menu item to specify the desired menu language.


- Touch the [right] arrow key in the entry, and then touch the [function enabled] icon of the desired language in the list of all available languages. The icon changes to [function activated]. The setting is immediately applied to the system.

### Show System information

Use this menu item to open the [System Information] window and load available updates. The following information is displayed:

- [Device id]
- [HMI Version]
- [HMI Update Version]: This entry changes to [Insert USB] stick as soon as you have inserted a USB stick with an HMI update into the USB slot on the bottom of the control panel.
- [Controller Version]
- [Controller Update Version]: This entry changes to [Insert USB stick] as soon as you have inserted a USB stick with a system update into the USB slot on the bottom of the control panel.
- Touch the [right] arrow key in the entry to open the [System Information] window.

When updating the HMI or system, do not remove the USB stick from the USB slot until the installation is complete. Note the progress indicator on the display.

When the installation is complete, a message appears at the bottom of the display. After each update, the device control restarts automatically.

## Datalog Settings

Use this menu to specify which system messages are to be written to the log file. All messages marked with the [function activated] icon are written to the log file. In the following example, its information texts, warning and error messages.

Datalog settings	
Log information	<input checked="" type="checkbox"/>
Log warnings	<input checked="" type="checkbox"/>
Log errors	<input checked="" type="checkbox"/>

## Show Datalog Series

Touch the arrow button [right] in the entry. The available measurement series are displayed. Measurement series are recorded using the Datalogging function, see Section "Data logging" on page 113.

## Error history

Touch the arrow button [right] in the entry to open the error history. All messages are displayed with a time stamp (example illustration).

- Touch the button [Remove individually] of an entry to remove it from the list.
- Touch the button [Remove all] at the bottom edge of the display to remove the whole list.

Error History	
	28-05-2015 14:38:17
	Bath Lift Sensor failure 28-05-2015 14:19:54
	The Bath Overheat Protection is switched on 28-05-2015 14:06:46
	Communication failure between Panel and Controller 28-05-2015 14:06:27
	Bath Lift Sensor failure 28-05-2015 14:05:07
	Communication failure between Panel and Controller 27-05-2015 15:17:58
	The Flask Support Rods are not in starting position

## Vacuum control mode

The device can be operated with manual (SETpressure) and automatic (AUTOaccurate) vacuum control.

For correct vacuum control, the vacuum valve, the vacuum sensor and the ventilation valve must be connected to the rear of the device (see the associated assembly instructions).

Regardless of the vacuum control mode, the vacuum pressure can be adjusted at any time during operation with the [Vacuum override] rotary control, see also section "Password / Screen lock" on page 97!

## SETpressure mode

In this operating mode, all process parameters are specified manually by the operator. To open the start window of this mode, touch the [Set Pressure] button on the start window:

### Index

<b>1</b>	SETpressure activated
<b>2</b>	Actual coolant temperature value from PT-100 sensor
<b>3</b>	[Save] button
<b>4</b>	Process parameter [Vacuum], actual/setpoint, start/stop button
<b>5</b>	Process parameter [Rotation], actual/setpoint display, start/stop button
<b>6</b>	Process parameter [Bath temperature], actual/setpoint, start/stop button
<b>7</b>	Process parameter [Steam temperature], actual value
<b>8</b>	[Stop all] button
<b>9</b>	[Warm-up] button
<b>10</b>	[Start all] button

- Touch a parameter's field to open the on-screen keyboard and set a desired value. The input values are accepted without a security prompt:
  - Vacuum setpoint (Adjustment Range: 1 – 1,400 mbar)
  - Hysteresis setpoint (1 – 50 mbar)
  - Rotation speed (6 – 160 rpm) (see also "Set rotation speed" on page 107)
  - Heating bath temperature water (20 – 100 °C) or oil (20 – 180 °C) (see also "Set the heating bath temperature" on page 108)

### **WARNING**

Incorrect/insufficient process settings can result in the risk of reaching unacceptable operating values: Risk of personal injury/property damage due to overheating, overpressure and/or escaping bath liquid/distillation material!

Only qualified personnel may adjust the process settings of the device.

When determining the pressure values, observe the chemical-physical properties of the respective material to be distilled.

Pay particular attention to the specific safety instructions for the processing of toxic substances!

A start/stop button pair is assigned to the parameters [Vacuum], [Rotation] and [Bath Temperature]: Touch these buttons to activate or deactivate the assigned process (build vacuum/maintain level, start/stop rotation, bath heating on/off) individually.

Touch the [Warm-up] button to activate the system warm-up function.

Touch the [Start All] button to start all processes at the same time with the displayed setpoints.

Touch the [Stop All] button to stop all processes at the same time.

The current settings can be stored and loaded at any time as Favorite in the system memory using the [Save] button (see section "Favorites and Profiles" on page 109).

## **AUTOaccurate mode**

In this operating mode, the vacuum pressure is regulated automatically and depending on the temperature. The rotation speed and heating bath temperature are manually set by the operator as in the SETpressure mode.

### **AUTOaccurate function**

The temperature-dependent vacuum control feature AUTOaccurate is based on two values:

- Inlet temperature cooling liquid: This value is measured with a temperature sensor (PT 100) at the chiller inlet on the rear of the device (see section "Mechanical design" on page 74)
- AUTOaccurate temperature: This temperature is measured at the AUTOaccurate sensor in the condenser

To use the function, a value  $\Delta T$  must be defined in the device control (parameter [Auto Accurate], see figure below). This value is added to the inlet temperature of the cooling liquid.

The sum  $\Delta T + \text{input temperature}$  of the cooling liquid serves as a reference value for the shutdown of the vacuum pump. This means that as soon as this temperature value is measured at the AUTOaccurate sensor in the condenser after the process has started, the evacuation stops.

As soon as the temperature at the AUTOaccurate sensor has subsequently decreased by a defined value again, the vacuum pump is switched on and the system will be evacuated again until the reference value at the sensor is reached again (with relatively increased vacuum pressure) and the vacuum pump switches off again.

In this way, a consistent performance is achieved in the processing of solvents or mixtures throughout the distillation process.

The described process steps are repeated until a defined maximum vacuum pressure is reached (input value [Vacuum & Steam temperature] → [end p]). The distillation process is then switched off.

**PRACTICAL TIP**

Before starting the process in AUTOaccurate mode, ensure that all set process temperatures are reached in a stable manner (observe display values!). Only then can the best possible and constant distillation results be achieved!

A warm-up time of 15 minutes is recommended as a guide. Use the device's warm-up function ([Warm Up] button in the footer of the start page).

- To open the start window of this mode, touch the [Auto Accurate] button on the start window:

**Index**

<b>1</b>	AUTOaccurate enabled
<b>2</b>	Actual coolant temperature value from PT-100 sensor
<b>3</b>	[Save] button
<b>4</b>	Process parameter [Max. vacuum pressure], actual value of the steam temperature
<b>5</b>	Process parameter [Rotation], actual/setpoint display, start/stop button
<b>6</b>	Process parameter [Bath temperature], actual/setpoint, start/stop button
<b>7</b>	Process parameter [ $\Delta T$ ], actual AUTOaccurate value
<b>8</b>	[Stop all] button
<b>9</b>	[Warm-up] button
<b>10</b>	[Start all] button

Touch a parameter's field to open the on-screen keyboard and set a desired value. The input values are accepted without a security prompt:

- Maximum vacuum pressure [end p] (setting range: 1 – 1,400 mbar)
- Rotation speed (6 – 160 rpm) (see also section "Set rotation speed" on page 107)
- Heating bath temperature water (20 – 100 °C) or oil (20 – 180 °C) (see also section "Set the heating bath temperature" on page 108)
- Input value  $\Delta T$  (0,1 – 30 °C)

### **WARNING**

Incorrect/insufficient process settings can result in the risk of reaching unacceptable operating values: Risk of personal injury/property damage due to overheating, overpressure and/or escaping bath liquid/distillation material!

Only qualified personnel may adjust the process settings of the device.

When determining the pressure values, observe the chemical-physical properties of the respective material to be distilled.

Pay particular attention to the specific safety instructions for the processing of toxic substances!

The hysteresis value for the AUTOaccurate operation cannot be adjusted. Factory setting: 0.5 °C.

The parameters [Rotation], [Bath Temperature] and [Auto Accurate] are assigned a start/stop pair of buttons. Touch these buttons to individually activate or deactivate the associated process (start/stop rotation, bath heating on/off, start/stop AUTOaccurate).

Touch the [Warm-up] button to activate the system warm-up function.

Touch the [Start All] button to start all processes at the same time with the displayed setpoints.

Touch the [Stop All] button to stop all processes at the same time.

The current settings can be stored and loaded at any time as Favorite in the system memory using the [Save] button (see section "Favorites and Profiles" on page 109).

### **Distillation under atmospheric pressure**

In certain cases, distillation under atmospheric pressure (low boiling) may be necessary. Refer to the specific safety and operating instructions in this section.

During atmospheric pressure distillation, only the flask rotation and the heating bath are started manually in the [SETpressure] operating mode after filling the evaporation flask.

The vacuum function must remain switched off, otherwise the vacuum pump will be switched on and vacuum will be built up in the system.

**CAUTION**

During distillation under atmospheric pressure, excess pressure may occur in the system. There is a risk that, in a high system pressure situation, gas or liquid will escape into the ambient air in the area of the seals and joints. Risk of personal injury/property damage

There is a risk of bursting under extreme pressure conditions!

Carefully compare the process parameters with the prevailing permissible ambient conditions and the chemical-physical properties of the respective material to be distilled (refer to data sheets)!

Pay particular attention to the specific safety instructions for the processing of toxic substances!

Only operate the device with the guard hood closed or under a closed laboratory fume hood.

Use appropriate personal protective equipment (PPE).

**Procedure**

- In the system settings, disable the [Leakage Pressure] stop criterion, see section "Security settings" on page 96. Touch the [function enabled] icon of the [Leakage pressure] entry. The icon changes to [function disabled].
- Open the start page again and switch to the SETpressure mode (see section "SETpressure mode" on page 102).
- Adjust the parameters [Rotation] and [Bath Temperature] taking into account the characteristics of the educt and the prevailing ambient conditions.
- Start the flask rotation and the heating bath using the associated start/stop buttons one after the other.

Safety Settings		
Overpressure	1400	mbar
Leakage Pressure	900	mbar <input checked="" type="checkbox"/>
Overheat protection	10	°C <input checked="" type="checkbox"/>

Do not start the vacuum function manually!

Do not start the distillation process with the [Start All] function, as this automatically activates the vacuum function.

Check the pressure in the evaporation flask at regular intervals (display value process parameter [Vacuum]).

## Set rotation speed

The procedure for adjusting the rotation speed is described in detail in sections "SETpressure mode" on page 102 and "AUTOaccurate mode" on page 103. Please also observe the information in this section!

---

### **WARNING**

When operating the device without a guard hood, there is a risk of contact with rotating parts. Loose clothing, jewelry and open hair can be pulled in!

- If possible, operate the device with the guard hood closed or under a closed laboratory fume hood.
- Never touch rotating parts with your hands. Do not bend over rotating parts.
- Wear tight-fitting work clothes, remove jewelry before starting work, tie long hair under a hair net or under a suitable head covering.

At high rotation speeds, the heating bath liquid can spray away from the evaporation flask by forming film and/or slop out of the heating bath. There is a risk of scalding and/or contamination.

- If possible, operate the device with the guard hood closed or under a closed laboratory fume hood.
- Adjust the rotation speed of the evaporation flask according to the recommendations in this section.
- Use appropriate personal protective equipment (PPE).
- Clean dirty surfaces immediately!

Even when the guard hood is closed, there is a risk that the heating bath liquid will spill out of the heating bath at high rotational speeds. There is a risk of scalding and/or contamination.

- Adjust the rotation speed of the evaporation flask according to the recommendations in this section.
- Use appropriate personal protective equipment (PPE).
- Clean dirty surfaces immediately!

## Recommended rotation speed

The speed of the evaporation flask affects the rate of distillation that can be achieved during the distillation process: The higher the rotation speed, the faster and more gentle the distillation material is processed.

Observe the following recommendations when using water as heating bath liquid:

Volume evaporation flask	Maximum recom- mended speed
20 liters	120 rpm
10 liters	140 rpm
6 liters	160 rpm

Due to the possible viscosity and volume range (depending on the oil type and the set heating bath temperature), no generally valid recommendation for the rotation speed of the evaporation flask can be given when using oil as heating bath liquid.

## Set the heating bath temperature

The procedure for adjusting the heating bath temperature is described in detail in sections "SETpressure mode" on page 102 and "AUTOaccurate mode" on page 103 . Please also observe the information in this section!

### CAUTION

The heating bath is designed for a maximum operating temperature of 210°C (using silicone oil). From a liquid temperature of 50 °C, there is a risk of injury in the event of contact.

At high rotational speeds, highly heated heating bath liquid can spray from the evaporation flask and/or swap out of the heating bath. There is a risk of scalding.


If possible, operate the device with the guard hood closed or under a closed laboratory fume hood.

Use appropriate personal protective equipment (PPE).

## Timer function

The timer function can be used to define a specific time period for a process sequence. At the same time, it is possible to specify which individual functions are stopped after a defined time has elapsed.

- Touch the [Timer...] button on the start window to open the [Set Timer] window:
- Touch the Days, Hours, Minutes and Seconds input fields to open the on-screen keyboard and enter the desired value (max. 99 days, 23 hours, 59 minutes, 59 seconds).
- In the selection [When Timer Ends:], specify which functions should be stopped after the timer expires (in the example above, all).
- Touch the [Close] button to dismiss all entries and return to the start window.

- Touch the [Start] button to start the timer and thus the process flow. The elapsed process time is visualized by means of a progress display on the start page:

When the timer expires, the message [Timer expired] appears. At the same time, the [Set Timer] window is opened again. Here you can read when the last process sequence and which individual functions were terminated.

The timer function can be stopped at any time with [Stop All]. In this case, the running process is canceled!

## Favorites and Profiles

The set process parameters can be stored in the system memory as [Favorite] (process parameter data set) in any operating mode. For the SETpressure operating mode, it is also possible to store certain process sequences as [Profile].

### Favorites mode

- Switch to the required operating mode and define the process parameters, see sections "SETpressure mode" on page 102 or "AUTOaccurate mode" on page 103.
- Touch the [Save] button:
- The [Save Favorite] dialog box appears. Touch the [Please enter name] input field and use the on-screen keyboard to set a name for the favorite.
- Confirm the entry with [Save] to save the new process parameter data set as a favorite.

The list of available favorites appears again.

or

- Finish the entry with [Confirm] to save the new process parameter data set as a favorite and load it immediately.
- Touch the [Delete] button to eliminate the opened favorite from the system memory.
- Touch the [Return] button to return to the list of available favorites without applying any changes.

### Create/edit favorites

The stored process parameter data records can be managed in the [Favorites] menu.

- Switch to the desired operating mode.
- From the start window, touch the menu button to open the main menu.
- Use [Favorites] to open the list of available favorites:
- Touch the [New] button to create a new record.
- or
- Touch the arrow button of an entry to open the corresponding record.
- Touch the [Edit] button in the footer to edit the new/existing favorite (see sections "SETpressure mode" on page 102 or "AUTOaccurate mode" on page 103).
- Load the record with [Apply].
- Touch the [Return] button to return to the list of available favorites without applying any changes.


## Load favorite

- Switch to the desired operating mode.
- From the start window, touch the menu button to open the main menu.
- Use [Favorites] to open the list of available favorites:
- Select the desired data record and confirm with [Apply].
- Then start the process with the loaded values with [Start All].

## Live view

- From the start window, touch the Show [Live View] button to open the process live view:

In the process live view, the parameters of the current process are visualized in real time:

The orange lines indicate the overall course.

The white colored area shows the previous progress (corresponding progress display).

The current process can be manually canceled at any time. To do this, touch the field with the name of the current favorite (at the top left of the start page) and confirm the security prompt.

## Create/edit process profiles

User-specific process sequences can be stored as profiles in the SETpressure operating mode:

- Switch to Set Pressure mode.
- From the start window, touch the menu button to open the main menu.
- Open the [Favorites] menu and switch to [Profiles] to open the list of available profiles:

### Edit an existing profile

- Touch the arrow button of an entry to open the corresponding record.
- Touch the [Edit] button in the footer to edit the profile (see also the following section "Create a new profile" on page 111).
- Confirm the changes with [Apply].
- Touch the [Return] button to return to the list of available profiles without applying any changes.


### Use existing profile as template

- Touch the arrow button of an entry to open the corresponding record.
- Touch the [Edit] button in the footer, and then touch [Copy] on the input screen to duplicate the profile.
- Use the on-screen keyboard to specify a label for the new profile and adjust the values as desired (see also the following section "Create a new profile" on page 111).
- Confirm the changes with [Apply].
- Touch the [Return] button to return to the list of available profiles without applying any changes.

### Create a new profile

- Touch the [New] button to create a new record.


- In the [Settings] input screen, first define a hysteresis value and the shutdown behavior for the individual functions (see section "Timer function" on page 108):

- Switch to the [Table] input screen:
- Define one process step per line here: Touch an input field (Time, Vacuum, Rotation, Bath Temperature) to open the on-screen keyboard.
- Set the desired value and confirm each entry with [OK].
- Insert new line: Highlight the line above which you want to insert a new line and touch the [New Line] button.
- Delete Row: Highlight the row you want to delete and touch the [Delete Row] button.

Set Pressure Profile abc			
Table	Settings	Graphs	
Time [hh:mm]	Vacuum [mbar]	Rotation [rpm]	Bath Temperature [°C]
0:00	800	0	20
0:01	800	60	60
0:02	800	80	80
0:03	1200	80	80
0:04	1200	120	90
0:05	1200	120	90

- After you have completed the parameterization, switch to the [Graphs] screen:
- Use the [ZOOM+/-] buttons to adjust the stretching of the time axes for the programmed ramp as desired.
- Confirm with [Save] to save the new profile in the system memory without loading it directly.
- Finish the entry with [Confirm] to save the new profile in the system memory and load it directly.
- Touch the [Return] button to return to the list of available profiles without applying any changes.

## Load process profile

- Switch to [Set Pressure] mode.
- From the start window, touch the menu button to open the main menu.
- Open the [Favorites] menu and switch to [Profiles] to open the list of available profiles:

Set Pressure Favorites	
★ Fixed Values	Profiles
abc 10-06-2015 11:46:18	▶
name 11-08-2015 11:16:39	▶

- Touch the arrow button of the desired profile and confirm with [Apply] to load the profile:
- Touch the [Return] button to return to the list of available profiles without applying any changes.

When the profile is completely processed, the message [Profile processed] appears. Here you can read when the profile was processed and which individual functions were terminated.

The processing of profiles can be canceled manually at any time with [Stop All].

## Data logging

With the datalogging function, the Hei-VAP Industrial process data can be recorded in the form of measurement series.

Data is recorded at individually defined intervals. With the shortest recording interval of one second, a maximum total recording time of approximately 27 hours can be achieved. If the Hei-VAP Industrial is operated in combination with an automatic module Hei-VOLUME Distimatic Pro, a recording interval of  $\geq 30$  seconds is recommended to ensure a recording period of at least three weeks.

When combined with an automatic module Hei-VOLUME Distimatic Pro, the Datalogging function of the Hei-VAP Industrial records the values of the integrated pressure sensor, including the overpressure values reached when draining the distillate out of the evaporator.

While draining the distillate by overpressure, which is controlled by the automatic module Hei-VOLUME Distimatic Pro, the process area (evaporation flask & condenser) is completely isolated by two closing valves and no medium is fed in or discharged.


During this period of up to 30 seconds, the vacuum values of the process area are not recorded!

After draining the distillate, the associated area is brought back to the process vacuum and the bypass is switched on again.

## Activating/deactivating data logging

The data logging function is manually activated/deactivated. The function status active/inactive is visualized using the appearance of the [Datalog...] button:


- Function inactive
- Function active

- Touch the [Datalog...] button to open the main window of the data logging function.

Depending on the operating state, one of the following status messages is displayed:

- [Data recording in progress]: Deactivate the datalogging function with [Stop]
- [Data recording stopped]: Activate the datalogging function with [Start].

The remaining storage time is displayed in days, hours, minutes, and seconds.

When all memory slots are occupied, the [Slots] button changes color. In this case, existing measurement series must be swapped out or deleted, see section "Deleting measurement series" on page 115.

- Touch the [Close] button to close the window.

## Setting the logging rate


The recording rate can be set in the main window:

- Touch the data entry fields [Data logging every XX minutes / seconds] one at a time and enter the desired values using the on-screen keyboard.

## Logging capacity

A total of 20 slots are available for recording a maximum of 100,000 measuring points. When the maximum recording capacity is reached, the following message appears:

- Confirm the message and proceed as described in the following section "Delete measurement series" to delete measurement series.

## Data list

The data list can be opened from the main window of the data logging function / button [Slots] or via the [Settings] / [Show Datalog Series] menu. The data list contains all available measurement series. Each entry consists of the time stamp and the number of measuring points recorded.

## Deleting measurement series

To delete an individual entry, touch the corresponding delete key (symbol **x**) and confirm the security prompt.

To delete all entries, touch the [Remove All] button in the footer and confirm the prompt.

To delete entries, the data logging function must be deactivated!

During the deletion process of all entries no individual deletion or saving processes are possible and the data list cannot be closed.

Datalog Series	
Start: 27-01-2021 8:50:25	<b>x</b>
End: 27-01-2021 8:51:29	
5 datapoints	
Start: 26-01-2021 16:29:33	<b>x</b>
End: 26-01-2021 16:42:06	
754 datapoints	
Start: 26-01-2021 15:52:58	<b>x</b>
End: 26-01-2021 16:42:19	
2960 datapoints	
Start: 26-01-2021 15:33:45	<b>x</b>
End: 26-01-2021 16:42:19	
4113 datapoints	

## Saving measurement series on a USB flash drive

Recorded measurement series can be saved on a USB stick. To do so, connect a VFAT formatted USB flash drive to the USB interface at the underside of the display and open the data list as described.

As soon as a suitable USB stick is connected, a [Save] button is assigned each list entry. At the same time, the [Save all] button appears in the footer of the main window.

Select a single entry in the list and touch the associated [Save] button to save that entry individually to the connected USB flash drive.

Alternatively, all list entries can be saved simultaneously.  
To do this, touch the [Save all] button in the footer:

During the saving process of all entries no individual saving processes are possible and the data list cannot be closed.

### Progress display when saving measurement series

When the saving process begins, the [Save] button is highlighted in orange. Once the process is completed successfully, the USB symbol will be replaced by a check mark.

The progress indicator is reset when you exit the mask.

## Data output

All data are saved in CSV and JSON format on the connected USB flash drive.

When opening the CSV file in Microsoft Excel®, the temperature values may be displayed in date format!

In this case, open the File / Options / Advanced menu and deselect the [Apply separator from operating system] option.

## Troubleshooting

System message	Possible remedy
Heating bath heating cannot be started	Water level too low, bath heating cannot heat up: Refill water bath.
Emergency off has been triggered	Press the emergency off button again, check the parameters and adjust if necessary, restart the process.
Rotation lock active	Rotary knob of the rotation lock is set to lock. Unlock.
Flask support not in start position	Support bars not in end position. Place support bars in end position.
Maximum pressure reached	In the vacuum control mode AUTOaccurate the value "end p" was reached, process was stopped. Distillation process finished!
Guard hood open	Guard hood was opened during a running process. Close the guard hood.
Communication between panel and controller interrupted	Internal COM error. Contact technical service!
Loading	Device operation while charging is in progress. No operator action until the message disappears.
Short circuit AUTOaccurate sensor	Contact technical service!
Overheat protection triggered	Check the level of the heating bath and the setting of the bath liquid selector on the rear of the device. If the heating bath heating unit is defective, contact authorized dealer or technical service.
Bath lift motor failure	Contact technical service!
Bath lift sensor failure	Contact technical service!
Bath mode (water/oil) cannot be changed during operation!	Bath liquid selector was actuated during a process. Move selector to original position.
Set bath temperature exceeded	Reduce the bath temperature and then increase gradually, clean the heating bath. If the heating bath heating unit is defective, contact authorized dealer or technical service.
Short circuit bath temperature sensor	Contact technical service!
Short circuit cooling water temperature sensor	Contact technical service!
Maximum permissible pressure exceeded	Decrease temperature and/or vacuum and increase gradually again, check setting in the Safety Settings menu.
Vacuum sensor broken	Replace vacuum sensor.
Short circuit vacuum sensor	Contact technical service!
Boiling temperature sensor not connected	Connect the boiling temperature sensor.
Short circuit boiling temperature sensor	Replace the boiling temperature sensor if defective.
AUTOaccurate sensor not connected	Contact technical service!
Bath heating temporarily suspended. Chiller temperature too high	Connect the AUTOaccurate sensor.
Bath temperature sensor broken	Replace the AUTOaccurate sensor if it is defective.
Coolant temperature sensor not connected	Reduce bath temperature, decrease hysteresis value, connect more powerful chiller.
Rotation lock activated during operation	Contact technical service!
	Connect the coolant temperature sensor.
	Replace coolant temperature sensor if defective.
	Unlock the rotation lock.

System message	Possible remedy
System leaking	Check the seals and connections, check the correct fit of the PTFE sealing, check the cleanliness and fit of the joints, grease if necessary, check the vacuum line and pump, check the setting in the Safety Settings menu.
Flask drive motor defective. Error Number:....	Contact technical service!
Actual vacuum value in the ventilated state differs from the air pressure.	Installation location is above normal zero, calibrate vacuum sensor.
	Flask support engaged. Slide outward.
	Guard hood open. Close.
Heating bath lift does not move up or down	Lift blocked. Check for blocking objects.
	Gap between housing and heating bath not free. Check for blocking objects.
	If the lift drive is defective, contact technical service.
	Flask support engaged. Slide outward.
No flask rotation	Rotation lock is in lock position. Unlock.
	Guard hood open. Close.
	If the flask drive is defective, contact technical service.
Ventilation valve does not switch	Connect the ventilation valve. If the valve is defective, contact authorized dealer or technical service.
Distillation prematurely stopped	Check end p value, check Set Time parameter.
Vacuum is not built up	Check PTFE seal, replace if necessary. Check sandwich seal, replace if necessary.

All error messages and warning messages are shown in plain text on the display. Follow the instructions on the display.

In case of recurring errors, please contact the responsible sales department or our technical service. Contact information see section "Contact details" on page 122.

## Technical specifications

### General device data

<b>Model</b>	<b>Rotary evaporator Hei-VAP Industrial</b>
Dimensions (W × H × D)	998 × 2,225 × 694 mm
Weight	approx. 120 kg, without glassware
Acoustic pressure	< 85 (dB(A)) (in accordance with IEC 61010)
Flask drive	EC motor
Speed range	6 – 160 rpm
Lift travel	180 mm
Lift speed	45 mm/s
Evaporation rates (L/h) ΔT* 40 °C (in continuous operation)	Toluene 26.0 Acetone 25.9 Ethanol 12.9 Water 4.7
Permissible leakage rate	1 mbar/min
Maximum cooling surface	1.2 m <sup>2</sup>
Control panel	detachable, 7 inch touchpanel, control knobs
Overtemperature protection	integrated, mechanical
Degree of protection, device	IP20
Control panel degree of protection	IP45

### Electrical data

Rated voltage	EU: 3/N/PE AC 230/400 V 50 Hz UNITED STATES: 2/PE AC 240 V 60
Connection	L1+L2+L3+N+PE or L+N+PE
Protection class	I (IEC 61140)
Oversupply category	II
Degree of pollution	2
Power input	4,600 W max.

### Heating bath, cooling

Heating bath	Inner Ø: 397 mm, outer Ø: 400 mm Material: V4A steel (1.4404) X2CrNiMo 17-12-2 Volume: 18 L H <sub>2</sub> O / 14 L oil
Heating power	4,000 W
Temperature range heating bath	20 – 100 C° H <sub>2</sub> O / 20 – 180 C° oil
Accuracy heating bath control	± 1 C°
Heating bath temperature control	Microprocessor
Required pump capacity	1,5 – 3 (m <sup>3</sup> /h)
Consumption cooling liquid	150 – 250 (L/h)

### Permissible ambient conditions

Storage conditions (recommendation)	5 °C – 31 °C at up to 80 % rel. humidity 32 °C – 40 °C at up to 50 % rel. humidity (decreasing linearly)
Operating temperature	5 °C – 31 °C at up to 80 % rel. humidity 32 °C – 40 °C at up to 50 % rel. humidity (decreasing linearly)
Installation altitude	up to 2,000 m asl

\*ΔT = difference between heating bath temperature and boiling temperature

### Scope of delivery

Item	Quantity	Product no.
Rotary evaporator Hei-VAP Industrial, complete*	1	518-13100-00
Warranty registration	1	01-006-002-78

\*The actual delivery depends on the respective device variant.

### Accessories

Detailed information on the available accessories for your device variant can be found on our website at [www.heidolph.com](http://www.heidolph.com).

### Device service

#### DANGER

Switch the device's main switch off and disconnect it from the power supply before carrying out maintenance work, cleaning, or repairs.

When cleaning, avoid the penetration of liquids.

Before replacing the fuses, switch off the device and disconnect the power supply cord.

Always replace the two device fuses in pairs with original manufacturer fuses. Further information on available accessories can be found on our website at [www.heidolph.com](http://www.heidolph.com)!

After fuse replacement, check the device for a safe condition according to IEC 61010-1.

### General cleaning instructions

#### CAUTION

Clean the device's surfaces with a soft, lint-free and only slightly moistened cloth.

Never use any aggressive or abrasive cleaning agents or aids.

### Cleaning and care instructions for heating baths

If water is used as a heating medium, it is generally recommended to change the bath liquid regularly.

In case of high salt load, use a mixture of 50 % tap water and 50 % deionized water to prevent calcification of the heating bath.

When using deionized or distilled water: Add 0.2 % borax ( $\text{Na}_2\text{B}_2\text{O}_7 \cdot 10 \text{H}_2\text{O}$ ) to deionized or distilled water.

If the heating bath container shows calcification, add a teaspoon of citric acid to the bath liquid before changing and heat the bath slightly until all the lime residue is dissolved or can be easily removed with a non-abrasive sponge. Then change the bath liquid.

Corroded areas can be carefully cleaned with a little stainless steel polish.

### Repairs - Return of equipment

Repairs to the device may only be carried out by authorized experts! Unauthorized repairs during the warranty period will result in the loss of the warranty claim. Regardless of the warranty claim, the owner is solely liable for damage caused by unauthorized repairs.

- In case of repair and before returning your device, contact our technical service at the following e-mail address:
  - [service@heidolph.de](mailto:service@heidolph.de).
- In your message, please provide us with the following information in addition to an error description:
  - Item number
  - Serial no.

The required data can be found on the rating plate of the device.

A service representative will contact you as soon as possible to agree on the next steps.

Include the completed certificate of decontamination with every device return, see "Certificate of decontamination" on page 123.

### Maintenance

There are no user-serviceable components in the unit housing. If necessary (in the event of abnormal operating behavior such as excessive noise or heat generation, for example), contact our technical service, see section "Contact details" on page 122.

- Check the caps and fasteners of all glass components for correct and tight fit at regular intervals and especially before each start-up of the device.
- Check the level sensor regularly for proper operation, contamination, and limescale build-up to prevent malfunctions.
- Check the PTFE seal of the EASY LOCK flask clamping system at regular intervals for heavy contamination and visible damage.
- Clean the PTFE seal of the EASY LOCK flask clamping system at regular intervals to prevent excessive deposits or encrustations.

---

#### CAUTION

At loose/leaky connections, ambient air is sucked in during the vacuum build-up and the required vacuum pressure cannot be reached in a stable manner!

In the event of excess pressure, there is a risk that liquid will escape from the system! The device and surrounding areas can be severely contaminated by leaking media!

Tighten loose caps/fasteners by hand or correct the seating of the retaining clips.

---

### Disposal

- When disposing of the device, observe the provisions of the WEEE Directive 2012/19/EU and its transposition into national law in the country of use.
  - When disposing of portable batteries, observe the provisions of the European Battery Directive 2013/56/EU and their implementation in national law in the country of use.
  - Check the device and all components for residues of substances that are hazardous to health, the environment and biohazardous before disposal.
  - Properly remove and dispose of residues of substances that are hazardous to health, the environment and biohazardous!
-

### **Warranty statement**

Heidolph Scientific Products GmbH provides a three-year warranty against material and manufacturing defects.

Glass and wear parts, transportation damage, and damage resulting from improper handling or non-intended use of the product are excluded from the warranty.

The warranty period for registered products begins on the date of purchase. Register the product with the enclosed warranty card or on our homepage [www.heidolph.com](http://www.heidolph.com).

For non-registered products, the warranty period begins with the date of the serial production (to be determined by the serial number).

In the event of material or manufacturing defects, the product will either be repaired or replaced free of charge within the warranty period.

### **Contact details**

#### **Heidolph Scientific Products GmbH**

Technical service

Walpersdorfer Str. 12

D-91126 Schwabach/Germany

Email: [service@heidolph.de](mailto:service@heidolph.de)

#### **Representations**

To find your local Heidolph distributor please visit [www.heidolph.com](http://www.heidolph.com)

### Certificate of decontamination

Enclose the certificate of decontamination, duly completed, with your device return. Submissions without a certificate of decontamination cannot be processed!

### CERTIFICATE OF DECONTAMINATION IN CASE OF RETURNS

Please fill in the required fields.

**Note: The sender must package the goods  
properly and appropriately for transport.**

Heidolph Scientific Products GmbH

Walpersdorfer Straße 12  
91126 Schwabach

Phone: +49 (0) 9122 9920-380  
**Fax: +49 (0) 9122 9920-19**  
E-Mail: service@heidolph.de

#### SENDER

Name \_\_\_\_\_

First name \_\_\_\_\_

Company/institution \_\_\_\_\_

Department \_\_\_\_\_

Address \_\_\_\_\_

Workgroup \_\_\_\_\_

ZC/City \_\_\_\_\_

Phone \_\_\_\_\_

Country \_\_\_\_\_

Email \_\_\_\_\_

#### DEVICE DETAILS

Article number \_\_\_\_\_

Serial no. \_\_\_\_\_

Ticket number \_\_\_\_\_

Reason for sending in \_\_\_\_\_

**Has the device been cleaned, decontaminated/disinfected?**

Yes

No

(Please mark as applicable)

If yes, which measures were carried out?

**Does this device pose a risk to people and/or the environment due  
to the processing of substances that are hazardous to health,  
the environment and/or are biohazardous?**

Yes

No

(Please mark as applicable)

If yes, with which substances did the device come into contact?

#### LEGALLY BINDING DECLARATION

The principal/consignor is aware that they are liable to the agent/consignee for losses or damage incurred due to incomplete and incorrect information.

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hereby declare, that the product (from serial number on 200067811) designated below is in compliance with the basic requirements of all applicable EU-directives stated below with regard to design, type of model sold and manufactured by us. This certificate will be invalid if the product is modified without the prior written consent and agreement of the manufacturer.

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A handwritten signature in blue ink, appearing to read 'J. K.' or 'J. H.'.

Andreas Hahn  
Managing Director

A handwritten signature in blue ink, appearing to read 'J. H.' or 'J. Z.'.

Jörg Ziel  
Quality Manager

**UK  
CA**

 **heidolph**  
research made easy

## Declaration of Conformity

In accordance with UK Government guidance

This declaration is issued under the sole responsibility of the manufacturer,  
Heidolph Scientific Products GmbH  
Walpersdorfer Straße 12  
91126 Schwabach / Germany

Product: Laboratory evaporator  
*Model:*

Hei-VAP Industrial xxx                            518-xxxxx-xx-x

Description:  
Rotary Evaporators with heating bath and human machine interface

The object of the declaration described above is in conformity with the relevant  
UK Statutory Instruments (and their amendments):

2008 No. 1597                                    *The Supply of Machinery (Safety) Regulations*  
2008 2016 No. 1091                                    *The Electromagnetic Compatibility Regulations 2016*  
2012 No. 3032    *The Restriction of the Use of Certain Hazardous Substances in  
Electrical and Electronic Equipment Regulations 2012*

and complies with the following technical standards :  
EN ISO 12100:2010, EN 61326-1:2013, EN 61010-1:2010, EN IEC 63000:2018

UK Authorised Representative (for authorities only):  
ProductIP ( UK ) Ltd.  
8. Northumberland Av.  
London WC2N 5BY

Signed for and on behalf of Heidolph Scientific Products GmbH  
Walpersdorfer Straße 12, 91126 Schwabach / Germany

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With this declaration, we confirm (according to current knowledge) that all sold laboratory devices and components by Heidolph Scientific Products GmbH fulfill the requirements of the EU directive 2011/65/EU (RoHS) and the extension 2015/863. All devices are compatible with the requirements of the RoHS Directive for the following materials:

Max. 0,01% of weight in homogenous material for cadmium and max. 0,1% of the weight in homogeneous material for Lead, Mercury, Hexavalent chromium, Polybrominated biphenyls (PBB), Polybrominated diphenyl ethers (PBDE), Bis(2-ethylhexyl) phthalate (DEHP), Butyl benzyl phthalate (BBP), Dibutyl phthalate (DBP), Diisobutyl phthalate (DIBP).

In the case of individual assemblies, maximum concentrations maybe exceeded within the permissible exceptions of the Directive.

Schwabach, 04.September.2024

Andreas Hahn  
Managing Director

Jörg Ziel  
Quality Manager

## China RoHS DECLARATION OF CONFORMITY

Heidolph Scientific Products GmbH has made reasonable efforts to ensure that hazardous materials and substances may not be used in its products.

In order to determine the concentration of hazardous substances in all homogeneous materials of the subassemblies, a "Product Conformity Assessment" (PCA) procedure was performed. As defined in GB/T 26572, the "Maximum Concentration Value" limits (MCV) apply to these restricted substances:

- Lead (Pb): 0.1%
- Mercury (Hg): 0.1%
- Cadmium (Cd): 0.01%
- Hexavalent chromium (Cr(VI)): 0.1%
- Polybrominated biphenyls (PBB): 0.1%
- Polybrominated diphenyl ether (PBDE): 0.1%

### Environmental Friendly Use Period (EFUP)

EFUP defines the period in years during which the hazardous substances contained in electrical and electronic products will not leak or mutate under normal operating conditions.

During normal use by the user such electrical and electronic products will not result in serious environmental pollution, cause serious bodily injury or damage to the user's assets.

The Environmental Friendly Use Period for Heidolph Scientific Products GmbH products is 25 years.

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This table is created according to SJ/T 11364-2014

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O: Indicates that the above mentioned hazardous substance contained in all homogeneous materials of the part is below the required limit as defined in GB/T 26572.

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X: Indicates that the above mentioned hazardous substance contained in at least one of the homogeneous materials of this part is above the required limit as defined in GB/T 26572.

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Apart from the disclosures in the above table, the subassemblies are not intentionally manufactured or formulated with lead (Pb), mercury (Hg), cadmium (Cd), hexavalent chromium (CrVI), polybrominated biphenyls (PBB), and polybrominated diphenyl ethers (PBDE).

Products manufactured by Heidolph Scientific Products GmbH may enter into further devices or can be used together with other appliances.

With these products and appliances in particular, Heidolph Scientific Products GmbH will not take responsibility for the EFUP of those products and appliances.

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Subject to change without notice. The printed version of this document is not regularly updated. The latest issue of this document can be found by visiting our homepage.

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MATERIAL CONTENT DECLARATION FOR Heidolph Scientific Products GmbH PRODUCTS
© Heidolph Cientific Products GmbH