

# PRECISE HEATING & STIRRING HEI-PLATE MIX 'N' HEAT MAGNETIC STIRRER

HIGH-RESISTANT HEATING PLATE

REFRACTORY HOUSING

800 W

MADE IN  
GERMANY

3-year warranty on  
all devices and an op-  
erational lifespan of  
more than 10 years

## Hei-PLATE Mix 'n' Heat Core

### Interface for more precision

For the EKT Hei-Con to measure and control the temperature precisely

### Robust and durable for all challenges in the lab

Kera-Disk® plate – scratch-resistant, refractory, chemical-resistant and easy to clean

Innovative body design for additional protection against burns, secure handling and reliable trapping of liquids from the device

Refractory die-cast aluminum housing

### Intuitive handling for a quick start

Standby button for switching the device easily on and off

Direct access to the main functions by two separate rotary controls

Heating button for the deliberate activation/deactivation of the heating function as well as for the visual display of the residual heat

Separate LED display of the heating



## ACCESSORIES RECOMMENDATION



Silicone protective cover



Electronic contact thermometer



Heating and bath attachments



High-performance condenser