

## 2540ET - 120V

### Technical datasheet



|                            |                                                       |
|----------------------------|-------------------------------------------------------|
| <b>Height</b>              | 8.3"                                                  |
| <b>Power</b>               | 120V, 50-60 Hz                                        |
| <b>Chamber Volume</b>      | 21L                                                   |
| <b>Chamber Dimensions</b>  | 10" x 18"                                             |
| <b>External Dimensions</b> | 18 x 17.7 x 23.8"                                     |
| <b>Wastewater Tank</b>     | 3.75L                                                 |
| <b>Connectivity</b>        | USB port (x2), Wifi compatibility, Ethernet port (x1) |
| <b>Number of Trays</b>     | 4                                                     |
| <b>Warranty</b>            | 3 Year                                                |
| <b>Certifications</b>      | UL, cUL, and ASME                                     |