

3850ELP - 230V

Technical datasheet



Height	12.1" (Max. Available)
Power	230V, 1 phase (60Hz) NEMA 6-15
Chamber Volume	65 L
Chamber Dimensions	15" x 19.7"
External Dimensions	28.3" x 21.3" x 30.1"
Utilities	No plumbing required
Autoclave Stand	023212609 (Included)
Autoclave Stand with Castors	023212601
Table - SS Top (36" x 48") Castors & Lower Shelf	023212602
Warranty	3 Year
Certifications	UL, cUL, and ASME