

## Hei-Dro Clean 165UC Washer

### Technical datasheet



|                                                    |                                             |
|----------------------------------------------------|---------------------------------------------|
| <b>Recommended Flow Pressure</b>                   | 44 psi (Min. 29 psi to Max. 73 psi)         |
| <b>Cold Water Temperature</b>                      | Min. 5°C - Max. 20°C                        |
| <b>Hot Water Temperature</b>                       | Min. 45°C - 70°C                            |
| <b>RO/DI Water Temperature</b>                     | Min. 5°C - Max. 20°C                        |
| <b>Connection Thread (with flask gasket)</b>       | G 3/4" (BSPP 3/4")                          |
| <b>Length of Water Inlet Hose</b>                  | 1.85 m/6.07'                                |
| <b>Max Drain Water Temperature</b>                 | 93°C/199°F                                  |
| <b>Max Pressure</b>                                | 2.9 psi                                     |
| <b>Max Drain Pump Head Height from Floor Level</b> | 1 m/3.28'                                   |
| <b>Chamber Volume</b>                              | 165 L                                       |
| <b>Chamber Dimensions</b>                          | (W x H x D) 21.7" x 19.7" x 23.6"           |
| <b>External Dimensions</b>                         | (W x H x D) 23.6" x 25.6" x 33.3"           |
| <b>Utility Requirements</b>                        | Cold Water, Hot Water, RO/DI Water<br>Drain |