

## Model 5075 ELVCP - 230V

### Technical datasheet



|                                                     |                                                         |
|-----------------------------------------------------|---------------------------------------------------------|
| <b>Height</b>                                       |                                                         |
| <b>Power</b>                                        | 230V, 3 phase (60Hz) NEMA L15-30                        |
| <b>Chamber Volume</b>                               | 160 L                                                   |
| <b>Chamber Dimensions</b>                           | 19.7" x 29.5"                                           |
| <b>External Dimensions</b>                          | 34.3" x 42.9" x 30.3"                                   |
| <b>Utilities</b>                                    | Treated Water, Tap Water, Compressed Air,<br>Floor Sink |
| <b>Wastewater Tank</b>                              | N/A                                                     |
| <b>Cooling Water Tank</b>                           | N/A                                                     |
| <b>Quiet Air Compressor for Models with Cooling</b> | 023250023                                               |