

System 27

Technical datasheet



Height	7 in/180 mm
Reaction Positions	24 (MTO format)
Stirring Speed	0 rpm to 1,100 rpm
Colume Range	0.5 ml to 7 ml
Temperature Range, Cooling	Ambient down to -7°C
Temperature Range, Heating	Ambient up to 150°C
Time to Max. Heat Up	30 min
Temperature Stablility	+/- 0.5°C
Base Diameter	6.9 in/176 mm
Top Diameter	7.5 in/190 mm
Well Diameter	0.49 in/12.5 mm
Operating Pressure	3 psi