

## Technical datasheet



|                                   |                      |
|-----------------------------------|----------------------|
| <b>Height</b>                     | 7 in/180 mm          |
| <b>Reaction Positions</b>         | 24 (MTO format)      |
| <b>Stirring Speed</b>             | 0 rpm to 1,100 rpm   |
| <b>Colume Range</b>               | 0.5 ml to 7 ml       |
| <b>Temperature Range, Cooling</b> | Ambient down to -7°C |
| <b>Temperature Range, Heating</b> | Ambient up to 150°C  |
| <b>Time to Max. Heat Up</b>       | 30 min               |
| <b>Temperature Stablility</b>     | +/- 0.5°C            |
| <b>Base Diameter</b>              | 6.9 in/176 mm        |
| <b>Top Diameter</b>               | 7.5 in/190 mm        |
| <b>Well Diameter</b>              | 0.49 in/12.5 mm      |
| <b>Operating Pressure</b>         | 3 psi                |