

Hei-SHAKE Rock Core 2000 (5°)

Technical datasheet



Weight	approx. 14 kg
Maximum permissible load	max. 15 kg
Drive	EC motor
Speed range	0 - 80 rpm
Acoustic pressure	< 85 dB (A), in accordance with IEC 61010
Rated voltage	1/N/PE 100-240V AC 50/60 Hz
Protection class	IP32 (acc. to EN 60529)
Power input	90 W
Movement	Rocking
Interfaces	1 × RS232 (Sub-D, 9 pole)
(W x H x D) with Platform 1000	n.a.
(W x H x D) with Platform 2000*1)	405 × 240 × 430 mm
Usable area (W x D) with Platform 1000	n.a.
Usable area (W x D) with Platform 2000	aprox. 390 × 340 mm
Tilt Angle	10°
Operating and Storage Temperature	5 °C – 31 °C at up to 80 % rel. humidity 32 °C – 40 °C at up to 50 % rel. humidity (decreasing linearly)
Installed Altitude	up to 2,000 m asl