

Technical datasheet



Hauteur	12.1" (Max. Available)
Puissance	230V, 3 phase (60Hz) NEMA L15-30
Chamber Volume	85 L
Chamber Dimensions	15" x 27.2"
External Dimensions	28.3" x 21.3" x 37"
Utilities	No plumbing required
Autoclave Stand	023212609 (Included)
Autoclave Stand with Castors	023212601
Table - SS Top (36" x 48") Castors & Lower Shelf	023212602
Warranty	3 Year
Certifications	UL, cUL, and ASME