

Hei-Dro Clean 165UC Washer

Technical datasheet



Recommended Flow Pressure	44 psi (Min. 29 psi to Max. 73 psi)
Cold Water Temperature	Min. 5°C - Max. 20°C
Hot Water Temperature	Min. 45°C - 70°C
RO/DI Water Temperature	Min. 5°C - Max. 20°C
Connection Thread (with flask gasket)	G 3/4" (BSPP 3/4")
Length of Water Inlet Hose	1.85 m/6.07'
Max Drain Water Temperature	93°C/199°F
Max Pressure	2.9 psi
Max Drain Pump Head Height from Floor Level	1 m/3.28'
Chamber Volume	165 L
Chamber Dimensions	(W x H x D) 21.7" x 19.7" x 23.6"
External Dimensions	(W x H x D) 23.6" x 25.6" x 33.3"
Utility Requirements	Cold Water, Hot Water, RO/DI Water Drain