

Duo™ II.I

Technisches Datenblatt



Auto Volumetric Dispense	0.1 - 25 L
Dimensions (Remote Dispenser)	640 x 250 x 255mm / 25.2" x 9.84" x 10"
Dispense Flow Rate	≥2 L/min
Water Quality (at 25° C)	18.2MΩ cm
Bacteria	<0.01 EU/ml
Make-Up Flow Rate (@15°C/59°F)	≥25 L/hr
Conductivity	<1400µS/cm
Collid Index	< 3 SDI
Dimensions (H x W x D)	500 x 490 x 290 mm / 19.7" x 19.3" x 11.5"
Drain Type	Sink/Tundish/Standard Up-Stand with Suitable Trap
Drain Flow	>3 L/min
DNASE	<5pg/mL