

Geno™ II

Technisches Datenblatt



Make-Up Flow Rate (@15°C/59°F)	≥25L/hr
Rejection Rate	≥98%
Rate from Tank	3 - 7 L/min
Typical TOC Reduction	96%
Recommended Daily Usage	160 L/day
Pressure	0.1 - 6 Bar/1.45 - 86 PSI
Conductivity	<1400µS/cm
Temperature	4 - 30°C/39.2 - 86°F
Hardness as CaCO3	<400 ppm
Free Chlorine	<2 ppm
Collid Index	<3 SDI
pH Value	4 - 10
Iron (Total Fe)	<0.1 ppm
Dimensions (H x W x D)	500 x 490 x 290 mm/19.7" x 19.3" x 11.5"
Installed Weight	23 kg/50.7 lbs
Power Supply	100 - 240 v ~; 50/60 Hz; 3.0A (max)
Power Supply Outbound	24 VDC, 5.0 A, 120 W (max)
Operating Temperature	4°C - 35°C max)/39.2°F - 95°F (max)
Drain Type	Sink/Tundish/Standard Up-Stand with Suitable Trap
Drain Flow	>3 L/min