

# Findenser™ Mini NS24 Cone, NS24 Socket

## Technical datasheet



|                                |                    |
|--------------------------------|--------------------|
| <b>Weight</b>                  | 670 g              |
| <b>Length × Diameter</b>       | 275 mm × 72 mm     |
| <b>Joints</b>                  | B14, B19, B24      |
| <b>Inner or wetted parts</b>   | Borosilicate glass |
| <b>Finned jacket</b>           | Anodised aluminium |
| <b>Plug &amp; Seal</b>         | Acetal & FKM       |
| <b>Thermofluid</b>             | Water              |
| <b>Max solvent bpt</b>         | 155 °C             |
| <b>Oven drying</b>             | Max Temp. 60 °C    |
| <b>Dishwasher</b>              | Max Temp. 50 °C    |
| <b>Autoclaving</b>             | Do Not Autoclave   |
| <b>Ambient operating temp.</b> | 0 °C to 24 °C      |